




# Theater Slim PCR Cycler

*User friendly, hot lid  
automatically heats up  
and it's portable*

## Features

- Power
  - 12V DC universal power supply or 12V DC battery pack
- Saves Space
  - Thin, only 4 cm
  - Very small footprint
  - Can be used in clean bench
- Touch Screen
- Saves Time
  - No waiting time
  - Finish setup in 30s

## Performance



PCR Product 500 bp

Amplifying Condition: 94°C 5 mins; 94°C 20 sec,  
56°C 30 sec, 72°C 1 min, 30 cycles; 72°C 10 mins



## Specifications

Capacity	0.2 ml*16
Temperature Range	8 - 99°C
Time Range	1s - 60 min
Display Precision	0.1°C
Temperature Precision	±0.2°C
Temperature Evenness (minimize edge effect)	±0.2°C
Number of Stored Program	15
Number of Cyclers	1 - 99
Control Interface	Touch Screen / Turning Knob
Dimensions	200*200*40 mm
Net Weight	2.1 kg
Max.Ramping Rate	3°C /s
Heated Lid Temperature	105°C (Can be switched off)
Operation Ambient Temperature	8 - 30°C
Power Supply	12V DC


## Applications

- Regular PCR
- Touch Down PCR
- Incubation
- RT-PCR
- Large Fragment PCR

## Set up your PCR reaction anywhere



- Clean Bench



- Electrophoresis Laboratory



- Bench

### Coyote Bioscience Co., Ltd.

2th Floor , Shangdixinilu#12 Haidian District, Beijing , China

Phone: +8610 64844237

[www.coyotebio.com/en](http://www.coyotebio.com/en)

### Coyote Bioscience USA Inc.

1506 Dell Avenue, Suite A, Campbell, CA 95008

Phone: +1 650-528-2199

[www.coyotebio-usa.com](http://www.coyotebio-usa.com)